Walgreens

PREMIUM ARM

Blood Pressure Monitor

- 5-Year Limited Warranty

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TABLE OF CONTENTS

INTRODUCTORY INFORMATION:

Important Product Notices and Safety Instructions	3
About Blood Pressure	5
Blood Pressure Standard.....	6
How This Blood Pressure Monitor Works	7

USAGE INFORMATION:

Name/Function of Each Part.....	8
Display Explanations	9
Installing Batteries	11
Using the AC Adapter.....	12
Date & Time Set Procedure.....	13
<i>Bluetooth</i>[®] Operation	15
Using Triple-Check (Multi-Read) Function	16
Using the UltraSoft[®] Arm Cuff	17
Check-with-Confidence[™]	21
Measurement Procedure.....	22
Transmit Readings To Your Mobile Device	26

Earn Balance Rewards Points.....	28
Risk Category Index	29
Irregular Heartbeat Detector	30
Excessive Body Motion Detector	31
Recalling Values from Memory.....	32
Clearing Values from Memory	35
Important Notes Regarding Your Blood Pressure Measurement.....	36
Care, Maintenance, & Cleaning	37
Potential for Electromagnetic Interference	39
Electromagnetic Compatibility (EMC).....	40
Troubleshooting.....	42
Specifications	43
<i>Bluetooth</i>[®] Specifications	45
5-Year Limited Warranty.....	46

IMPORTANT PRODUCT NOTICES AND SAFETY INSTRUCTIONS

When using your blood pressure monitor, basic precautions should always be followed. Please read and follow all instructions and warnings before using this product. Save these instructions for future reference.

- This is a home health care product only and it is not intended to serve as a substitute for the advice of a physician or medical professional.
- This device uses the oscillometric method to measure systolic and diastolic blood pressure, as well as heart rate.
- **DO NOT** use this device for diagnosis or treatment of any health problem or disease. Measurement results are for reference only. Consult a health care professional for interpretation of pressure measurements. Contact your physician if you have or suspect any medical problem. Do not change your medications without the advice of your physician or health care professional.
- Proper cuff size is critical for accurate measurements. Follow the instructions in this manual and printed on the arm cuff to ensure the appropriate size cuff is being used.
- For those who have undergone mastectomy surgery (especially those who have had lymph nodes removed), it is recommended to take a measurement on the unaffected side.
- When using medical electronic equipment(s) on the same limb, pressurization of the cuff may cause temporary malfunction to the other devices.
- Wearing the arm cuff properly is critical for accurate measurements. If the circumference of the upper arm differs largely from the circumference of the forearm, an excessive gap may occur between the arm and the arm cuff, which could result in an inaccurate measurement. Please consult your healthcare professional if you have questions on how to properly wear the arm cuff.
- This product is not suitable for people with arrhythmias, those undergoing intravenous injections on any limb, or those currently on dialysis treatment. This device may have difficulty determining the proper blood pressure for users with irregular heartbeat, diabetes, poor circulation of blood, kidney problems, or for users who have suffered from a stroke.

- This monitor is clinically validated for use during pregnancy and preeclampsia. If you experience unusual readings, even after the recommended rest period, please consult your medical professional.
- Excessive use may result in blood flow interference, which is likely to cause uncomfortable sensations, such as partial subcutaneous hemorrhage, or temporary numbness to your arm. In general, these symptoms should not last long. However, if symptoms persist, please seek advice from a medical professional.
- The pulse display is not suitable for checking the frequency of heart pacemakers.
- Electromagnetic interference: The device contains sensitive electronic components. Avoid strong electrical or electromagnetic fields in the direct vicinity of the device (ex. mobile telephones, microwave ovens). These may lead to temporary impairment of measurement accuracy.
- Consider the electromagnetic compatibility of the device (ex. power disturbance, radio frequency interference, etc.) Please use this device in a home health care environment only.
- Use blood pressure monitor only for its intended use.
- **DO NOT** wrap the cuff around body parts other than your arm.
- **DO NOT** touch, move, or shake the device during a measurement, as this could cause an inaccurate blood pressure reading.
- **DO NOT** talk or move your arm or hand muscles during a measurement, as this could cause an inaccurate blood pressure reading.
- **DO NOT** use this device on infants, children, or those who cannot express their own intention.
- Not for use by or on persons under the age of 18.
- Use only 1.5V “AA” alkaline batteries for power supply.
- The patient is an intended operator.
- The applied part is the cuff.
- Blood pressure measurements determined with this device are equivalent to those obtained by a trained observer using the cuff/stethoscope auscultatory method, within the accuracy limits prescribed by the American National Standard, manual, electronic, or automated sphygmomanometers. The sphygmomanometer was clinically investigated according to the requirement of ISO 81060-2: 2013.

ABOUT BLOOD PRESSURE

Blood pressure is the pressure exerted on the artery walls while blood flows through the arteries. The pressure measured when the heart contracts and sends blood out of the heart is systolic (highest) blood pressure. The pressure measured when the heart dilates with blood flowing back into the heart is called diastolic (lowest) blood pressure. Both the systolic and diastolic pressure are necessary for a physician to evaluate the status of a patient's blood pressure.

Why Measure Your Blood Pressure?

Among today's various health problems, those associated with high blood pressure are very common. High blood pressure dangerously correlates with cardiovascular diseases. Therefore, blood pressure monitoring is important for identifying those at risk.

Why Do My Readings Vary?

Blood pressure is a body parameter that is subject to normal variations throughout the day. A single reading that is different from yours or your doctor's readings is not necessarily inaccurate. The average of several readings, taken under similar conditions, using the same arm is preferred for accurate blood pressure readings.

Why Are My Readings Different Than Those Taken at My Doctor's Office?

Many experience a phenomenon called "White Coat Hypertension" when measured by a doctor. White Coat Hypertension refers to blood pressure that rises above its usual level when measured in a clinical setting, such as a doctor's office.

BLOOD PRESSURE STANDARD

This table contains defined levels for hypertension that are publicly available from the American Heart Association® (AHA 2017). Users can compare their own blood pressure readings against these defined levels to determine if they may potentially be at an increased risk.

This table is applicable to most adults aged 18 and older.

Blood Pressure Category	Systolic mmHg (upper number)		Diastolic mmHg (lower number)	LED Indicator Color
Normal	<120	and	<80	Green
Elevated	120–129	and	<80	Yellow
High Blood Pressure (hypertension) Stage 1	130–139	or	80–89	Red
High Blood Pressure (hypertension) Stage 2	140–180	or	90–120	
Hypertensive Crisis (consult your doctor immediately)	>180	and/or	>120	

*Source: AHA 2017

Blood pressure tends to go up and down, even in people who normally don't have high readings. If your numbers stay above the normal range most of the time, you may be at increased risk and should consult your physician.

Although one can easily find where their own blood pressure readings fall on this table, this monitor comes equipped with a Risk Category Index that automatically compares each reading to the defined levels and provides a helpful cue if your reading falls into one of the stages that could potentially indicate increased risk. See page 29 for more information on this feature.

Please note that cues provided by this monitor are only intended to assist you in using this table. The table and cues are only provided for convenience to help you understand your noninvasive blood pressure reading as it relates to the American Heart Association information. They are not a substitute for a medical examination by your physician. It is important for you to consult with your physician regularly. Your physician will tell you your normal blood pressure range as well as the point at which you may actually be considered to be at risk.

HOW THIS BLOOD PRESSURE MONITOR WORKS

This monitor uses innovative GentleInflate™ technology to detect your blood pressure. With one touch of a button, the monitor automatically takes the reading during inflation, then deflates, delivering a quick, comfortable and quiet measurement every time.

Please note that any muscle movement during inflation will cause measurement error. When measurement is complete, the monitor will display your systolic pressure, diastolic pressure, and pulse readings.

The monitor automatically finds where your measurement results fall on the American Heart Association® (AHA 2017) table and provides a cue if your reading falls into one of the stages that could potentially indicate increased risk. Please refer to page 29 for more information on this feature.

The appearance of the "♥ icon indicates that a pulse irregularity consistent with an irregular heartbeat was detected during measurement. Refer to page 30 for more information on the Irregular Heartbeat Detector.

NAME/FUNCTION OF EACH PART

NOTE:


The **“MEM”** button is also used to set the date and time.

Other Accessories:

4 “AA” 1.5V alkaline batteries (included)


DISPLAY EXPLANATIONS



Display Symbols:

	User 1: Appears when the monitor is operated by User 1.
	User 2: Appears when the monitor is operated by User 2.
	Low Battery Symbol: Appears when batteries should be replaced, or +/- polarities have been positioned incorrectly.
	Pulse Symbol: Shows the heart rate per minute.
	Irregular Heartbeat Detector: See page 30 for more information.
	Risk Category Indicator: The bar points out the specific risk category that your reading falls into. The color of the bars will vary, depending on which category your reading falls into. The bar(s) will be either green (normal), yellow (elevated), or red (high).


	<p>AM/PM Average: Indicates the reading being displayed is an average from the last 3 morning or last 3 evening measurements.</p>
	<p>Data Transmitting: Appears when data is transferring to your mobile device.</p>
	<p>Triple-Check (Multi-Read) Function Result: Indicates which measurement is being taken, or which measurement is being viewed from a Triple-Check (Multi-Read) average reading.</p>
	<p>Triple-Check (Multi-Read) Function Symbol: Appears when Triple-Check (Multi-Read) function is turned ON.</p>
	<p>Triple-Check Average Symbol: Displayed when viewing a Triple-Check average.</p>
	<p>Excessive Body Motion Detector: Appears when excessive body movement, especially of the arm the cuff is worn, is detected during the measurement. See page 31 for more information. NOTE: The measured blood pressure reading may not be accurate if the icon is displayed. Properly wrapping the arm cuff is important. If the circumference of the upper arm differs largely from the circumference of the forearm, an excessive gap may occur between the arm and the arm cuff, and result in the false detection of excessive body motion.</p>
	<p>Check-with-Confidence™: Appears upon first use and then every 30th use, during the self-check function. See page 21 for more information.</p>

If **ERROR** and any of the following letters and numbers appear in the area that systolic pressure should be displayed, an error has occurred with your reading. See Troubleshooting section of this manual for more information.

	<p>Measurement Error: Wrap the cuff correctly and keep arm steady during measurement. Measure again.</p>
	<p>Air Circuit Abnormality: Check cuff connection. Measure again.</p>
	<p>Pressure Exceeding 300 mmHg: Turn the unit off to clear, then measure again.</p>
	<p>Error Determining Measurement Data: Wrap the cuff properly, keep steady, and measure again.</p>
	<p>Data Transmission Error: Monitor cannot connect to the mobile device to transmit data. Make sure <i>Bluetooth</i>® is ON.</p>
	<p>System Error: Turn off monitor and measure again. If "EP" error still appears on the display, please call a Consumer Relations representative.</p>

INSTALLING BATTERIES

1. Press down on latch and lift the battery cover to open the battery compartment.
2. Install or replace 4 “AA” batteries in the battery compartment according to the indications inside the compartment.
3. Close the battery cover by pushing in the top end of the battery door.

Replace the batteries if:

- The low battery symbol appears on the display.
- When any button is pressed and nothing is displayed on the screen.

NOTE:


- Date and time will need to be reset if batteries are removed or replaced.
- Replace all batteries at one time (as simultaneous set). Do not mix old and new batteries.
- Use only 1.5V “AA” alkaline batteries. Do not mix alkaline, standard (carbon-zinc), or rechargeable (Ni-Cad, Ni-MH, etc.) batteries.
- When installing batteries, observe proper +/- polarities. Incorrect installation of batteries may cause damage to the unit.
- Remove batteries when unit is not in use for extended periods of time.
- When the batteries are removed, the measurement values stored in memory are retained. However, the date and time must be reset.
- Clean contacts on battery and in battery compartment with a soft, dry cloth each time you install batteries.
- Recycle or dispose of properly in accordance with all local, state, province, and country regulations.
- Keep batteries away from children as they could pose a choking hazard.
- **Batteries are hazardous waste. Do not dispose of them together with the household garbage.**
- **DO NOT dispose of batteries in fire. Batteries may explode or leak.**

USING THE AC ADAPTER

1. Connect the AC adapter with the AC adapter jack as shown below.
2. Please use only the AC adapter included with this monitor.

NOTE:

- Date and time will need to be reset if AC adapter is unplugged and unit is without batteries.
- When the AC adapter is your main power supply, make sure the adapter plug can be easily removed from the unit.
- When using with the adapter, remove the batteries when not in use for an extended period of time. Leaving the batteries in the compartment for an extended period of time may cause leakage, which may lead to damage to the unit.
- No batteries are needed when operating with an AC adapter.


Connect AC
Adapter Here

AC adapter

Input: 100-240V~0.18A-0.20A 50/60Hz

Output: 6V DC 1A 6W Max

DATE & TIME SET PROCEDURE

Date and time can be set by two methods, either sync automatically using your *Bluetooth*® mobile device, or set manually using the “**Date & Time Set**” ⌚ and “**MEM**” buttons.


If setting the date and time by syncing with your mobile device, it is important that this is done prior to taking any measurements to ensure the date and time are accurate.

To set using your device with Bluetooth wireless technology:

1. Before taking the first measurement, open the **BP Health App** on your device. If you have not already installed the **BP Health App** onto your device, it is available to download on the App StoreSM and on Google PlayTM. Make sure you have the app downloaded and open on your mobile device before trying to sync the date and time.

2. Press the **Bluetooth**  button on the monitor to enable the connection.

3. Once the *Bluetooth* connection is established, the date and time will automatically update on the blood pressure monitor.

NOTE:


- Once the date and time have been successfully synced, future readings will automatically have the correct date and time.
- If you are having trouble automatically syncing the date and time, make sure the *Bluetooth* is ON on your device and ON on the blood pressure monitor.

To Set Manually:

1. To set the date and time, press the **“Date & Time Set”** button .

2. The display will show a blinking number showing the MONTH. Change the MONTH by pressing the **“MEM”** button. Each press will increase the number by one in a cycling manner. Press the  button again to confirm the entry and the screen will show a blinking number representing the DAY.
3. Change the DAY, YEAR, HOUR, and MINUTE as described in step 2 above, using the **“MEM”** button to change the numbers and the  button to confirm the entries.

Repeat process to set DAY, YEAR, HOUR, and MINUTE.


IMPORTANT: IF YOU ARE USING THIS BLOOD PRESSURE MONITOR FOR THE FIRST TIME, PLEASE REMOVE THE PROTECTIVE FILM FROM THE DISPLAY SCREEN.

BLUETOOTH® OPERATION

This monitor has *Bluetooth* function ON for your convenience. This will allow your readings to automatically transmit to the **BP Health App**. By turning this function OFF, the measurements cannot be transmitted.

Turning *Bluetooth* Function OFF:


While the screen is off (nothing is shown on the screen), press and hold the **Bluetooth**  button or 3 seconds to turn *Bluetooth* OFF.

- When *Bluetooth* is OFF, measurements cannot be transmitted to the **BP Health App**.

Turning *Bluetooth* Function ON:

While the screen is off (nothing is shown on the screen), press and hold the **Bluetooth**  button for 3 seconds to turn *Bluetooth* ON.


- When *Bluetooth* is ON, readings can automatically and manually be transmitted to the app.
- The **BP Health App** must be open on your mobile device in order to establish a connection with the monitor.

Bluetooth QDID Profile Information:

QDID: 45008, 51582, & 52727

USING TRIPLE-CHECK (MULTI-READ) FUNCTION

The Triple-Check (Multi-Read) Function automatically takes and averages 3 readings in a row, with 1 minute rest intervals in between each measurement. To deactivate this feature and take only a single reading, slide the switch on the side to the OFF position.


USING THE ULTRASOFT® ARM CUFF

Very Important: Proper cuff size is critical to accurate measurement.

This monitor comes with one UltraSoft Upper Arm Cuff that fits arm circumferences from 9"-17" (23cm-43cm).


If the arrow falls outside of the solid-colored line, you will need a larger cuff. A larger cuff, designed for arms from 17"-22" (43cm-56cm), is available for this monitor. To purchase, please contact HoMedics Consumer Relations at 1-800-466-3342.


NOTE:

- If an oversized cuff is used, the blood pressure measurement tends to be lower than the actual blood pressure.
- If an undersized cuff is used, the blood pressure measurement tends to be higher than the actual blood pressure.

1. If the cuff is not assembled, pass the end of the cuff furthest away from the tubing through the metal D-ring in order to form a loop. The smooth side without the felt material should be on the inside of the cuff loop.

2. Plug the cuff tube into the left side of the unit.

3. Open the arm cuff so that your arm may be placed through it.
4. Remove tight-fitting clothing from your upper left arm so that the cuff can be wrapped around your exposed arm.




NOTE: It is recommended to wrap the cuff on a bare arm or over thin clothing to ensure accuracy. Thick clothing or a rolled up sleeve will cause inaccurate blood pressure measurements.

5. Position cuff edge 0.8-1.2 inches (2-3 cm) above elbow.


NOTE: If using the larger cuff designed for arms 17"-22" (43cm-56cm), position cuff above inner elbow crease.

6. Center tube over middle of arm.

7. Pull the end of the cuff so that it tightens evenly around your arm. Press the hook and loop material together to secure, ensuring the cuff does not slip during measurement. Allow room for 2 fingers to fit between the cuff and your arm.


8. Lay your arm on a table (palm upward) so the cuff is at the same height as your heart. Make sure the tube is not kinked.

NOTE:

- This device should not be used when your arm is wounded or injured.
- Thick clothing or a rolled up sleeve will cause inaccurate blood pressure measurements. It is recommended the cuff be wrapped on a bare arm or over thin clothing.
- If it is not possible to fit the cuff on your left arm, it can also be placed on your right arm. However, all measurements should be made using the same arm.
- To use on the right arm, you must position the artery symbol “ Φ ” over the main artery. Locate the main artery by pressing with 2 fingers approximately 1" (2 cm) above the bend of your elbow on the inside of your right arm. Identify where the pulse can be felt the strongest. This is your main artery.

CHECK-WITH-CONFIDENCE™

This monitor includes a built-in Check-with-Confidence feature provides accuracy assurance function, which enables the device to automatically perform a self-check. It provides a convenient way to check and ensure the monitor is working properly.

The Check-with-Confidence default setting is “ON.” To turn this feature “OFF”, press the **“START/STOP”** button for 3 seconds.

Check-with-Confidence
feature ON

Check-with-Confidence
feature OFF

On first use, then on every 30th use, the cuff will partially inflate and then deflate. The monitor will display “OK”, confirming the monitor is working properly.

If the monitor is not working properly, “Err” will be displayed on the screen. Please call Consumer Relations at 1-800-466-3342 for further assistance.

NOTE:

If the **“MEM”** and **“START/STOP”** buttons are pressed at the same time, erasing the values in memory, the monitor will perform the self-check feature on the next measurement.

MEASUREMENT PROCEDURE

IMPORTANT NOTES:

- Blood pressure changes with every heartbeat and is in constant fluctuation throughout the day.
- Blood pressure measurement can be affected by the position of the user, his or her physiologic condition, and other factors.

Before Measurement:

To help ensure an accurate reading, follow these instructions before taking a measurement:

- Wait 1 hour after exercising, bathing, eating, drinking beverages with alcohol or caffeine, or smoking to measure blood pressure.
- Sit quietly and rest for 15 minutes.
- Stress raises blood pressure. Avoid taking measurements during stressful times.
- Take your reading in a comfortable environment as measurements can be affected by hot or cold temperatures. Take your blood pressure at normal body temperature.

During Measurement:

- DO NOT talk or move your arm or hand muscles.
- DO NOT cross your legs. Sit with feet flat on the floor.
- DO NOT touch cuff or monitor during measurement.

Using Single Measurement Mode:

1. Please make sure the Triple-Check switch, located on the side of the unit, is in the OFF position.


2. Press the “1Δ/Δ2” button to choose User 1 or User 2.

After the user number is selected, press the “START/STOP” button to confirm the user.

3. With the cuff wrapped around your arm, press the “START/STOP” button. Do not inflate the cuff unless it is wrapped around your arm. All digits will light up to check the display functions. The checking procedure will be completed after about 3 seconds.

NOTE: The Check-with-Confidence™ icon  will stay on the screen if this is the 1st or 30th use, and the feature is “ON”. The monitor will then perform a self-check to ensure it is working properly.

4. After all symbols disappear, the display will show “00”. The monitor is ready to measure and will automatically inflate the cuff to begin measurement.
5. When the measurement is completed, the cuff will deflate entirely. Systolic pressure, diastolic pressure, and pulse will be shown simultaneously on the LCD screen. **The measurement is then automatically stored into memory.**

NOTE:

- This monitor automatically returns to standby mode approximately 1 minute after last operation. You may also press the “**START/STOP**” button to return to standby mode.
- To interrupt the measurement, you may press the “**START/STOP**” button (recommended) or “**MEM**” buttons. The cuff will deflate immediately after a button is pressed.

Using Triple-Check Measurement Mode:

1. Please make sure the Triple-Check switch, located on the side of the unit, is in the ON position.

2. Press the “**18/12**” button. After the user number is selected, press the “**START/STOP**” button to confirm the chosen user.
3. With the cuff wrapped around your arm, press the “**START/STOP**” button. Do not inflate the cuff unless it is wrapped around your arm. All digits will light up to check the display functions. The checking procedure will be completed after about 3 seconds.


NOTE: The Check-with-Confidence™ icon ✓ will stay on the screen if the feature is on.

4. After all symbols disappear, the display will show “00”. The monitor is ready to measure and will automatically inflate the cuff to begin measurement.
5. When the measurement is completed, the cuff will deflate entirely. Systolic pressure, diastolic pressure, and pulse will be shown simultaneously on the LCD screen. **The measurement is then automatically stored into memory.**


NOTE:

- If an error code is displayed during measurement, the monitor will repeat the measurement up to 3 times during the consecutive measurement process.
 - Remain still until all 3 measurements are completed.
6. At the conclusion of the measurement, the Triple-Check average will be displayed.

To review your individual results that make the Triple-Check average, press the **“Date & Time Set”** ⌚ button.

press the
⌚ button

TRANSMIT READINGS TO YOUR MOBILE DEVICE

This monitor syncs your blood pressure readings to the **BP Health App** that is free to download on the App StoreSM and on Google PlayTM. Make sure you have the **BP Health App** downloaded, and open on your mobile device before trying to transmit your blood pressure measurements.


Before attempting to sync the blood pressure monitor with your mobile device, make sure Bluetooth[®] is turned ON, on both your mobile device and the monitor.

Automatically Transmit Readings:

After a measurement is taken, the *Bluetooth* icon  will appear on the screen as the monitor automatically transmits your blood pressure readings to the App.

Manually Transmit Readings:

Press the *Bluetooth*  button. The readings will automatically be transmitted to the app. The Bluetooth icon will appear on the screen as the monitor transmits your blood pressure readings to the app.

If the transmission(s) were successful, the *Bluetooth* icon  will be displayed on the screen.

If the transmission(s) were unsuccessful, E4 and the **ERROR** symbol will be displayed on the screen.

NOTE:

- The measurement is stored in the monitor's memory even if the reading is not transmitted to your mobile device.
- The **BP Health App** must be open on your mobile device in order to transmit your measurements.
- Only new readings will be accepted by the app.
- The blood pressure monitor can only be assigned to one Walgreens user account.
- To ensure readings transmitted to the app have the correct date and time, it is important that the correct date and time is set on the blood pressure monitor before taking measurements.
- Measurements transmitted to the app cannot be edited.
- This monitor is subject to and complies with electromagnetic compatibility (EMC) standard of IEC 60601-1-2, EN 301 489-1, EN 301 489-17, EN 300 328 and U.S. federal guidelines, Part 15 of the FCC (Federal Communications Commission) rules for devices with RF capability. These guidelines help ensure that your monitor will not affect the operation of other nearby devices. Additionally, other devices should not affect the use of your monitor.
- Other wireless devices that are in use nearby, such as a cell or mobile phone, or a wireless network, may prevent or delay the transmission of data from your monitor to paired *Bluetooth*[®]-enabled device. Moving away from the source of the interference or turning off these devices may resolve the problem.
- Make sure the monitor and paired *Bluetooth*-enabled device are within an acceptable distance (no more than 30 feet) away from each other.
- Make sure to select the correct user on the monitor before your blood pressure measurement begins.
- *Bluetooth* data transmission is not available while taking a measurement.

EARN BALANCE REWARDS POINTS*

Earn Balance Rewards points for taking blood pressure readings when you sync this device to your smartphone and use the **BP Health App**.

*Balance Rewards

- Limitations apply to earning points.
- Information provided to Walgreens for the Balance Rewards program is subject to the terms and conditions of Balance Rewards and is protected in accordance with the Online Privacy and Security Policy.
- See full terms and details at [BalanceRewards.com/healthychoices](https://www.balance-rewards.com/healthychoices).

RISK CATEGORY INDEX

This monitor comes equipped with a Risk Category Indicator that automatically compares each reading to defined levels established by the American Heart Association® (AHA 2017) as described earlier in this manual, and provides a helpful cue if your reading falls into one of the stages that could potentially indicate increased risk. Please note that cues provided by this monitor are only intended to assist you in using this table.

The table is only provided for your convenience to help you understand your noninvasive blood pressure readings as it relates to the AHA (2017) information. They are not a substitute for a medical examination by your physician. It is important for you to consult with your physician regularly. Your physician will tell you your normal blood pressure range as well as the point at which you may actually be considered to be at risk.

Blood Pressure Category	Systolic mmHg (upper number)		Diastolic mmHg (lower number)	LED Indicator Color
Normal	< 120	and	< 80	Green
Elevated	120–129	and	< 80	Yellow
High Blood Pressure (hypertension) Stage 1	130–139	or	80–89	Red
High Blood Pressure (hypertension) Stage 2	140–180	or	90–120	
Hypertensive Crisis (consult your doctor immediately)	> 180	and/or	> 120	

*Source: AHA 2017

IRREGULAR HEARTBEAT DETECTOR


The appearance of the "♥️" icon indicates that a pulse irregularity consistent with an irregular heartbeat was detected during measurement. Usually this is not a cause for concern. However, if the symbol appears often, we recommend you seek medical advice. Please note that the device does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage.

Movement, shaking or talking during the measurement can result in pulse irregularities that may cause the appearance of this icon. Therefore, it is of great importance to not move or talk during measurement.

To determine the presence of an irregular heartbeat, the average of the heartbeat intervals is calculated with the first 3 normal effective heartbeat values. It is important to note that the average is not a strict mathematical averaging of all recorded intervals. At least 3 beats with 25% or greater difference from the average heartbeat interval will generate the "♥️" icon on the screen.

Important Information:


This blood pressure monitor is not designed for use by people with arrhythmias nor for diagnosing or treating an arrhythmia problem. As a safeguard, we recommend that if you have arrhythmias such as atrial or ventricular premature beats and atrial fibrillation or any other special conditions, you should check with your physician before using your blood pressure monitor.

Irregular Heartbeat Detector

EXCESSIVE BODY MOTION DETECTOR

The (⚡) icon is displayed if body movement is detected during measurement, especially the movement of the arm the blood pressure monitor is worn on.

Excessive Body Motion Detector



NOTE:

- The measured blood pressure reading may not be accurate if the icon is displayed.
- If you have purchased and are using the XL cuff for arm sizes 17"-22", a gap between the cuff and the arm could cause this icon to display, please adjust cuff as directed and measure again.

RECALLING VALUES FROM MEMORY

This monitor can be used by 2 individuals. Each user can store up to 120 measurements.

This monitor features an advanced memory mode to provide you with a variety of options for reviewing your measurement history, including AM/PM averaging. Reviewing your morning (AM) and nighttime (PM) measurements can provide important information about how your blood pressure changes throughout the day.

	Displays average of last 3 readings.
	Displays average of last 3 morning readings. (AM is defined as 4:00 AM – 12:00 PM).
	Displays average of last 3 nighttime readings. (PM is defined as 6:00 PM – 2:00 AM).
	Indicates reading is an average of a Triple-Check measurement. Press the “Date & Time Set” ⌚ button to review individual readings that make up the Triple-Check average.

NOTE:


A Triple-Check average and the 3 measurements that make this average, count as 4 readings when stored in memory.

To use the memory function:


1. Press the “1Δ/Δ2” button to select User 1 or User 2.

2. Press the “MEM” button to access the memory.


3. The monitor will first display the calculated average applied to the last 3 memories.

4. Press the “MEM” button again and the monitor will display the average of the last 3 AM memories.

NOTE: If there are fewer than 3 AM measurements in memory, this average will not be displayed.


5. Press the “MEM” button again and the monitor will display the average of the last 3 PM memories.

NOTE: If there are fewer than 3 PM measurements in memory, this average will not be displayed.


6. Every new press of the “MEM” button will recall a previous reading. The latest reading will be recalled first.

NOTE:


- If the last 3 readings were a Triple-Check average, then the Triple-Check average will be displayed.
- Memory averaging function will only average individual readings.

NOTE: While reviewing the readings, the Triple-Check average icon may appear on the screen. Press the “Date & Time Set” ⌚ button to review individual readings that make up the Triple-Check average.


CLEARING VALUES FROM MEMORY


1. Press “1Δ/Δ2” button to select User 1 or User 2.

2. Press the “MEM” button to enter into memory recall mode.

3. Press and hold the “MEM” and the “START/STOP” buttons at the same time and the data for the selected user will be erased automatically.

NOTE:

- Once deleted, your readings cannot be restored.
- Once deleted, the device will automatically perform a self-check on the next measurement.

IMPORTANT NOTES REGARDING YOUR BLOOD PRESSURE MEASUREMENT

- It is suggested that you take your measurements at the same time each day and use the same arm for consistency.
- Users should wait a minimum of 5 minutes before taking additional measurements. More time may be necessary depending upon your physiology.
- The measurement results that users receive are for reference only. If users have any blood pressure concerns, please consult a physician.
- Once inflation reaches 300 mmHg, the unit will deflate automatically for safety reasons.
- This product is not suitable for people with arrhythmias.
- This device may have difficulty determining the proper blood pressure for users with irregular heartbeat, diabetes, poor circulation of blood, kidney problems, or for users who have suffered a stroke.

CARE, MAINTENANCE, & CLEANING

- Clean the blood pressure monitor body and cuff carefully with a slightly damp, soft cloth. Do not press. Do not wash cuff or use chemical cleaner on it. Never use thinner, alcohol, or petrol (gasoline) as cleaner.
- Leaky batteries can damage the unit. Remove the batteries when the unit will not be used for a long time.
- Follow local ordinances and recycling instructions regarding disposal or recycling of the device and device components, including batteries.
- If the unit is stored near freezing, allow it to acclimate to room temperature before use.
- This blood pressure monitor is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems with this device, please contact HoMedics Consumer Relations (contact information can be found on the Warranty page).
- **DO NOT** immerse the unit in water as this will result in damage to the unit.
- **DO NOT** subject the monitor or cuff to extreme temperatures, humidity, moisture, or direct sunlight. Protect from dust.
- **DO NOT** fold the cuff and tubing tightly.
- **DO NOT** disassemble the monitor or cuff. If in need of repair, refer to the warranty section of this manual.
- **DO NOT** subject the monitor to extreme shock (do not drop on floor).
- **DO NOT** inflate the cuff unless wrapped around arm.
- **DO NOT** wrap the cuff around body parts other than your arm.
- **DO NOT** drop or insert any object into any opening or hose.

- To avoid accidental strangulation, keep this product away from children. Do not drape tube around neck.
- This monitor may not meet its performance specifications if stored or used outside of these temperature and humidity ranges:

Storage/Transportation Environment

Temperature:
-13°F~158°F (-25°C~70°C)

Humidity:
less than 93% RH

Operating Environment

Temperature:
41°F~104°F (5°C ~40°C)

Humidity:
15% ~ 93% RH

POTENTIAL FOR ELECTROMAGNETIC INTERFERENCE

To avoid inaccurate results caused by electromagnetic interference between electrical and electronic equipment, do not use the device near a cell phone or microwave oven. For most wireless communication devices, it is recommended to maintain a distance of 10.8 feet (3.3m) in order to avoid electromagnetic interference.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following 2 conditions: (1) This device may not cause harmful interference, and (2) This device must accept any interference received, including interference that may cause undesired operation.

The equipment complies with FCC RF exposure limits set forth for an uncontrolled environment. The equipment must not be co-located or operated in conjunction with any other antenna or transmitter.

FEDERAL COMMUNICATIONS COMMISSION COMPLIANCE STATEMENT

Changes or modifications to this equipment not expressly approved by the manufacturer could void the user's authority to operate the equipment.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.
- Any changes or modifications not expressly approved by the grantee of this device could void the user's authority to operate the equipment.

ELECTROMAGNETIC COMPATIBILITY (EMC)

• Guidance and manufacturer's declaration – electromagnetic emissions

The device is intended for use in the electromagnetic environments listed below, and should only be used in such environments:

Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	RF energy is used only to maintain device's operation. Therefore, its RF emissions are so low that it's not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments, and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

• Guidance and manufacturer's declaration – electromagnetic immunity

The device is intended for use in the electromagnetic environments listed below, and should only be used in such environments:

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment–guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact discharge ± 15 kV air discharge	± 8 kV contact discharge ± 15 kV air discharge	In air discharge testing, the climatic conditions shall be within the following ranges: Temperature: 15° – 35° , Humidity: 30% – 60%.
Power frequency (50 or 60 Hz) magnetic field IEC 61000-4-8	30 A/m 50 or 60 Hz	3 A/m 50 or 60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 2 kV Power lines	± 2 kV Power lines	Mains power quality should be that of a typical commercial or hospital environment.
Interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315° 0% UT; 1 cycles 70% UT; 25/30 cycles 0% UT; 250/300 cycles	0% UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315° 0% UT; 1 cycles 70% UT; 25/30 cycles 0% UT; 250/300 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.

• **Recommended separation distances between portable and mobile RF communication equipment and the device.**

The device is intended for use in an electromagnetic environment where radiated RF disturbances are under control. User can help prevent electromagnetic interference by keeping the device at a minimum distance from portable and mobile RF communications equipment (transmitters). The table below details the maximum output power of transmitter:

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.7 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

• **Guidance and manufacturer's declaration – electromagnetic immunity**

The device is intended for use in the electromagnetic environments listed below, and should only be used in such environments:

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3V rms At 0.15-80 MHz 6V rms At ISM & Radio Amateur Freq.	3V rms At 0.15-80 MHz 6V rms At ISM & Radio Amateur Freq.	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated RF IEC 61000-4-3 (Proximity fields from RF wireless communications equipment IEC 61000-4-3)	10 V/m at 80-2700 MHz AM Modulation And 9-28V/m at 385-6000MHz, Pulse Mode and other Modulation. The system shall be tested as specified in IEC60601-1-2 Table 9 for proximity fields from RF wireless communications equipment using the test methods specified in IEC 61000-4-3	10 V/m at 80-2700 MHz AM Modulation And 9-28V/m at 385-6000MHz, Pulse Mode and other Modulation. The system shall be tested as specified in IEC60601-1-2 Table 9 for proximity fields from RF wireless communications equipment using the test methods specified in IEC 61000-4-3	Recommended separation distance $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.7 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.

b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

TROUBLESHOOTING

If any abnormality arises during use, please check the following points:

SYMPTOMS	POSSIBLE CAUSES	CORRECTION
Unit does not turn on when the "START/STOP" button is pushed.	Batteries have run down.	Replace them with 4 new "AA" alkaline batteries.
	Battery polarities have been positioned incorrectly.	Re-insert the batteries in the correct positions.
EE measurement error symbol shown on display or the blood pressure value is displayed excessively low (or high).	Cuff has been placed incorrectly.	Wrap the cuff properly so that it is positioned correctly. Measure again.
E1 error symbol shown on display.	Air circuit abnormality. Cuff tube may not be plugged into monitor correctly.	Check cuff connection. Measure again.
E2 error symbol shown on display.	Inflation pressure exceeding 300 mmHg.	Switch the unit off, then measure again.
E3 error symbol shown on display.	Error determining measurement data.	Turn off monitor and measure again. If 'EP' error still appears on the display, please call Consumer Relations.
	If the circumference of the upper arm differs largely from the circumference of the forearm, an excessive gap may occur between the arm and the arm cuff.	Ensure cuff is wrapped correctly. Please consult your healthcare professional if you have questions on how to properly wear the arm cuff.
E4 data transmission error symbol shown on the display.	Bluetooth function is turned OFF, on your mobile device.	Turn Bluetooth ON, on your mobile device, and try again.
	Mobile device does not pair with the blood pressure monitor.	Make sure Bluetooth is ON, on the mobile device and the blood pressure monitor, and try again. Requires compatible mobile device. This blood pressure monitor is compatible with mobile devices that are running iOS 5.0 or newer, or Android 4.3 and newer.
	The App on the mobile device is not open.	Make sure the app is open on the mobile device, and try again.
	The blood pressure monitor and mobile device are out of transmitting range.	Make sure the mobile device and blood pressure monitor are within the acceptable range of 32 feet (10 meters).
	Unexpected loss of electrical/mechanical integrity.	Make sure your last reading is stored in memory and the App is open and try again. See page 32. Remove batteries, re-insert, and try again. Return the device to your local distributor or importer.
EP error symbol shown on display	System error.	Measure again.
 Excessive body motion detector symbol shown on display.	Excessive body movement, especially of the arm the cuff is worn, occurred during the measurement. Ex: talking, moving, or shaking of the arm with the cuff on while measurement.	Measure again. Keep arm steady during measurement.
	Cuff may not be wrapped correctly, or the circumference of the upper arm differs largely from the circumference of the forearm, an excessive gap may occur between the arm and the arm cuff.	Ensure arm cuff is wrapped correctly. See "Using the UltraSoft® Arm Cuff" on page 17 for more information. Please consult your healthcare professional if you have questions on how to properly wear the arm cuff.

NOTE: If the unit still does not work, contact HoMedics Consumer Relations. Under no circumstance should you disassemble or attempt to repair the unit by yourself. Contact information for HoMedics Consumer Relations can be found on the Warranty page.

SPECIFICATIONS

Power Source:	Four 1.5V DC "AA" alkaline batteries or 6V DC AC adapter
Measurement Method:	Oscillometric
Rated Range of Cuff Pressure	0~300 mmHg
Rated Range of Determination	40~280 mmHg
Measurement Range of Heart Rate	40~199 beats/minute
Accuracy:	Pressure: ± 3 mmHg Pulse: $\pm 5\%$ of reading
Pressure Sensor:	Semi-conductor
Inflation:	Pump driven
Deflation:	Automatic pressure release valve
Memory Capacity:	240 memory total for 2 users
Operation Environment:	Temperature: 41°F~104°F (5°C ~40°C) Humidity: 15% ~ 93% RH Atmospheric pressure: 700 hPa~1060 hPa
Storage/Transportation Environment:	Temperature: -13°F~158°F (-25°C~70°C) Humidity: Less Than 93% RH
Weight:	0.67 lbs (303g) (without batteries)
Arm Circumference:	One-Size UltraSoft® Cuff—Fits Most Arms: 9"–17" (23cm–43cm) XL Cuff (not included): 17"–22" (43cm–56cm)
Dimensions:	4.45"(L) x 5.24"(W) x 2.20"(H) 113mm(L) x 133mm(W) x 56mm(H)
Battery Life:	Approx. 200 measurements
Product Life:	5 years (4 measurements per day)
Accessories:	(4) "AA" alkaline batteries, 6V DC AC adapter, instruction manual, 1 arm cuff with tube
Auto Shutoff:	Without any operation for 1 minute, device automatically shuts off
RF Type:	Bluetooth® 4.0 BLE
System Requirement of the Bluetooth-enabled Device:	Bluetooth 4.0 for Android 4.3 or above Bluetooth 4.0 for iOS 7.0 or above

	<p>Internally powered equipment Type BF Equipment Not suitable for use in presence of flammable anesthetic mixture with air, oxygen or nitrous oxide. Continuous operation with short-time loading.</p>
	<p>Follow instructions for use.</p>
	<p>To avoid inaccurate results caused by electromagnetic interference. NOTE: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the device. Otherwise, degradation of the performance of this equipment could result.</p>
	<p>Waste of electrical and electronic equipment (WEEE)</p>
	<p>Serial number SN YYMMXXXXXX</p>
<p>IP22</p>	<p>Ingress Protection Rating First characteristic numeral – Degree of protection against access to hazardous parts and against foreign objects. N1=2 (Protected against solid foreign objects of 12.5 mm \varnothing and greater). Second characteristic numeral – Degree of protection against ingress of water. N2=2 (Protected against vertically falling water drops when ENCLOSURE tilted up to 15°).</p>
	<p>Non-ionizing electromagnetic radiation.</p>

NOTE: These specifications are subject to change without notice. If the device needs to be checked for calibration, please call Consumer Relations.

BLUETOOTH[®] SPECIFICATIONS

RF Type	<i>Bluetooth</i> 4.0 BLE
RF Modulation	GFSK
Equivalent Isotropic Radiated Power	-4.18dBm
Data Throughput	0.2Mbps
Expected Delay (Latency Range) in Wireless (RF) Communication	The latency time is less than 0.3ms from sender to receiver.
Integrity	Channel Quality-Driven Data Rate (CQDDR) technology increases the effective data rate and integrity in noisy environments.
Security	AES-128 and application layer user defined
Wireless Operation Distance	Class 2 (Maximum: 10 meter)
RF Frequency / Need for Spectrum Management	2402 - 2480 MHz (allowing for guard bands)
Maximum Limitation	Unlimited
Maximum Permitted Power	0.47mW
Proximity of Other In-band Transmitters Used in Vicinity	up to 40 bands (2 MHz spacing; centered from 2402 to 2480 MHz)
Wireless Communication Profile	GATT – Client and Server
Wireless Coexistence	Support for 802.11 Coexistence
System requirement of the Bluetooth-enabled device	<i>Bluetooth</i> 4.0 for Android 4.3 or above <i>Bluetooth</i> 4.0 for iOS 7.0 or above

Made exclusively for

Walgreens

**For service or repair,
do not return this unit
to the retailer. Contact
HoMedics Consumer
Relations at:**

**Email:
cservice@
homedics.com**

**Phone:
1-800-466-3342
Business Hours:
8:30am-7pm ET
Monday-Friday**

**Distributed by
HOMEDICS®**

**HoMedics, LLC
3000 N Pontiac Trail
Commerce Township, MI
48390**

Printed in China

LIMITED 5-YEAR WARRANTY

HoMedics sells its products with the intent that they are free of defects in manufacture and workmanship for a period of 5 years from the date of original purchase, except as noted below. HoMedics warrants that its products will be free of defects in material and workmanship under normal use and service. This blood pressure monitor meets the simulated measurement cycles test requirement per EN1060-3, part 8.10. This warranty extends only to consumers and does not extend to Retailers.

To obtain warranty service on your HoMedics product, contact a Consumer Relations Representative by telephone at 1-800-466-3342 for assistance. Please make sure to have the model number of the product available.

HoMedics does not authorize anyone, including, but not limited to, Retailers, the subsequent consumer purchaser of the product from a Retailer or remote purchasers, to obligate HoMedics in any way beyond the terms set forth herein. This warranty does not cover damage caused by misuse or abuse; accident; the attachment of any unauthorized accessory; alteration to the product; improper installation; unauthorized repairs or modifications; improper use of electrical/power supply; loss of power; dropped product; malfunction or damage of an operating part from failure to provide manufacturer's recommended maintenance; transportation damage; theft; neglect; vandalism; or environmental conditions; loss of use during the period the product is at a repair facility or otherwise awaiting parts or repair, or any other conditions whatsoever that are beyond the control of HoMedics. This warranty is effective only if the product is purchased and operated in the country in which the product is purchased. A product that requires modifications or adoption to enable it to operate in any other country than the country for which it was designed, manufactured, approved and/or authorized, or repair of products damaged by these modifications is not covered under this warranty.

THE WARRANTY PROVIDED HEREIN SHALL BE THE SOLE AND EXCLUSIVE WARRANTY. THERE SHALL BE NO OTHER WARRANTIES EXPRESS OR IMPLIED INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS OR ANY OTHER OBLIGATION ON THE PART OF THE COMPANY WITH RESPECT TO PRODUCTS COVERED BY THIS WARRANTY. HOMEDICS SHALL HAVE NO LIABILITY FOR ANY INCIDENTAL, CONSEQUENTIAL, OR SPECIAL DAMAGES. IN NO EVENT SHALL THIS WARRANTY REQUIRE MORE THAN THE REPAIR OR REPLACEMENT OF ANY PART OR PARTS WHICH ARE FOUND TO BE DEFECTIVE WITHIN THE EFFECTIVE PERIOD OF THE WARRANTY. NO REFUNDS WILL BE GIVEN. IF REPLACEMENT PARTS FOR DEFECTIVE MATERIALS ARE NOT AVAILABLE, HOMEDICS RESERVES THE RIGHT TO MAKE PRODUCT SUBSTITUTIONS IN LIEU OF REPAIR OR REPLACEMENT.

This warranty does not extend to the purchase of opened, used, repaired, repackaged and/or resealed products, including but not limited to sale of such products on internet auction sites and/or sales of such products by surplus or bulk resellers. Any and all warranties or guarantees shall immediately cease and terminate as to any products or parts thereof which are repaired, replaced, altered, or modified, without the prior express and written consent of HoMedics.

This warranty provides you with specific legal rights. You may have additional rights which may vary from state to state. Because of individual state regulations, some of the above limitations and exclusions may not apply to you.

©2019 HoMedics, LLC. All rights reserved.

HoMedics and UltraSoft are registered trademarks of HoMedics, LLC.

GentleInflate is a trademark of HoMedics, LLC.

American Heart Association is a registered trademark of the American Heart Association.

The Bluetooth word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Walgreens is under license.

App Store is a service mark of Apple, Inc.

The Apple logo is trademark of Apple, Inc., registered in the U.S. and other countries.

IOS is a registered trademark of Cisco in the U.S. and other countries and is used under license by Apple Inc.

Android and Google Play are trademarks of Google, Inc.

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©2019 HoMedics, LLC. Todos los derechos reservados.
HoMedics y UltraSoft son marcas registradas de HoMedics, LLC.
GentleInflate es una marca comercial de HoMedics, LLC.
App Store es una marca de servicio de Apple Inc.
El logotipo de Apple es una marca comercial de Apple, Inc., registrada en los EE. UU. y en otros países.
Android y Google Play son marcas registradas de Google, Inc.