Walgreens

DELUXE WRIST

Blood Pressure Monitor

- 5-Year Limited Warranty

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TABLE OF CONTENTS

INTRODUCTORY INFORMATION:

Important Product Notices and Safety Instructions	3
About Blood Pressure	5
Blood Pressure Standard.....	6
How This Blood Pressure Monitor Works.....	7

USAGE INFORMATION:

Name/Function of Each Part	8
Display Explanations	9
Installing Batteries	12
Date & Time Set Procedure.....	14
Applying the Cuff.....	15
Correct Measurement Posture	16
Wrist Positioning Guide	17
Wrist Positioning Guide Settings	18

Measurement Procedure.....	19
Risk Category Index.....	21
Irregular Heartbeat Detector	22
Excessive Body Motion Detector	23
Recalling Values from Memory.....	23
Clearing Values from Memory.....	25
Important Notes Regarding Your Blood Pressure Measurement.....	26
Care, Maintenance, & Cleaning	27
Potential for Electromagnetic Interference.....	29
Electromagnetic Compatibility (EMC).....	30
Troubleshooting.....	34
Specifications	36
5-Year Limited Warranty	38

IMPORTANT PRODUCT NOTICES AND SAFETY INSTRUCTIONS

When using your blood pressure monitor, basic precautions should always be followed. Please read and follow all instructions and warnings before using this product. Save these instructions for future reference.

- This device uses the oscillometric method to measure systolic and diastolic blood pressure, as well as heart rate.
- Please note that this is a home healthcare product only and it is not intended to serve as a substitute for the advice of a physician or medical professional.
- Not for use by or on persons under the age of 18.
- This device is only intended for use by those individuals with a wrist circumference ranging from 5.3"-8.4" (13.5cm-21.5cm).
- **DO NOT** use this device for diagnosis or treatment of any health problem or disease. Measurement results are for reference only. Consult a healthcare professional for interpretation of blood pressure measurements. Contact your physician if you have or suspect any medical problem. Do not change your medications without the advice of your physician or healthcare professional.
- This product is not suitable for people with arrhythmias.
- This device may have difficulty determining the proper blood pressure for pregnant women and for individuals with irregular heartbeat, diabetes, poor circulation of blood, kidney problems or for users who have suffered from a stroke.
- Not suitable for people undergoing intravenous injection on any limb, women who may be pregnant and those with preeclampsia.
- Those who have had mastectomy surgery (especially when lymph nodes have been removed), it is recommended to take measurements on the unaffected side.
- Excessive use may result in blood flow interference, which is likely to cause uncomfortable sensations, such as partial subcutaneous

hemorrhage, or temporary numbness to your wrist. In general, these symptoms should not last long. However, if symptoms persist, please seek advice from a medical professional.

- When used along with other medical electronic equipment on the same arm, pressurization of the cuff may cause the other devices to temporarily malfunction.
- Consult your physician before measuring blood pressure at the wrist if you have any of the following conditions: Severe hypertension, severe diabetes, severe arteriosclerosis, severe kidney disease, or vascular disease that may compromise circulation.
- The pulse display is not suitable for checking the frequency of heart pacemakers.
- Not suitable for use on patients with implanted electrical devices, such as a cardiac pacemakers or defibrillators.
- Electromagnetic interference: The device contains sensitive electronic components. Avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g. mobile telephones, microwave ovens). These may lead to temporary impairment of measurement accuracy.
- Consider the electromagnetic compatibility of the device (ex. power disturbance, radio frequency interference, etc.)
- Please use this device in a home health care environment only.
- Use blood pressure monitor only for its intended use.
- **DO NOT** wrap the cuff around body parts other than your wrist.
- **DO NOT** use this device on infants, children, or those who cannot express their own intention.
- Use only 1.5V “AAA” alkaline batteries for power supply.

Blood pressure measurements determined with this device are equivalent to those obtained by a trained observer using the cuff/stethoscope auscultatory method, within the accuracy limits prescribed by the American National Standard, manual, electronic, or automated sphygmomanometers. The sphygmomanometer was clinically investigated according to the requirement of ISO 81060-2:2013.

ABOUT BLOOD PRESSURE

What is Blood Pressure?

Blood pressure is the pressure exerted on the artery walls while blood flows through the arteries. The pressure measured when the heart contracts and sends blood out of the heart is systolic (highest) blood pressure. The pressure measured when the heart dilates with blood flowing back into the heart is called diastolic (lowest) blood pressure.

Why Measure Your Blood Pressure?

Among today's various health problems, those associated with high blood pressure are very common. High blood pressure dangerously correlates with cardiovascular diseases. Therefore, blood pressure monitoring is important for identifying those at risk.

Why Do My Readings Vary?

Blood pressure is a body parameter that is subject to normal variations throughout the day. A single reading that is different from yours or your doctor's readings are not necessarily inaccurate. The average of several readings, taken under similar conditions, using the same arm is preferred for accurate blood pressure readings.

Why Are My Readings Different Than Those Taken at My Doctor's Office?

Many experience a phenomenon called "White Coat Hypertension" when measured by a doctor. White Coat Hypertension refers to blood pressure that rises above its usual level when measured in a clinical setting, such as a doctor's office.

BLOOD PRESSURE STANDARD

This table contains defined levels for hypertension that are publicly available from the American Heart Association® (AHA 2017). Users can compare their own blood pressure readings against these defined levels to determine if they may potentially be at an increased risk.

This table is applicable to most adults aged 18 and older.

Blood Pressure Category	Systolic mmHg (upper number)		Diastolic mmHg (lower number)	Indicator Color
Normal	<120	and	<80	Green
Elevated	120–129	and	<80	Yellow
High Blood Pressure (hypertension) Stage 1	130–139	or	80–89	Red
High Blood Pressure (hypertension) Stage 2	140–180	or	90–120	
Hypertensive Crisis (consult your doctor immediately)	>180	and/or	>120	

*Source: AHA 2017

Blood pressure tends to go up and down, even in people who normally don't have high readings. If your numbers stay above the normal range most of the time, you may be at increased risk and should consult your physician.

Although one can easily find where their own blood pressure readings fall on this table, this monitor comes equipped with a Risk Category Index that automatically compares each reading to the defined levels and provides a helpful cue if your reading falls into one of the stages that could potentially indicate increased risk. See page 21 for more information on this feature.

Please note that cues provided by this monitor are only intended to assist you in using the table on the previous page. The table and cues are only provided for convenience to help you understand your noninvasive blood pressure reading as it relates to the American Heart Association® information. They are not a substitute for a medical examination by your physician. It is important for you to consult with your physician regularly. Your physician will tell you your normal blood pressure range as well as the point at which you may actually be considered to be at risk.

HOW THIS BLOOD PRESSURE MONITOR WORKS

This monitor uses innovative GentleInflate™ technology to detect your blood pressure. With one touch of a button, the monitor automatically takes the reading during inflation, then deflates, delivering a quick, comfortable and quiet measurement every time.

To help guide your wrist into the ideal position, this blood pressure monitor comes equipped with an optional wrist positioning guide to help you take an accurate blood pressure reading. Please see page 17 for more information on this feature.

Please note that any muscle movement during inflation will cause a measurement error. When measurement is complete, the monitor will display your systolic pressure, diastolic pressure, and pulse readings.


The monitor automatically finds where your measurement results fall on the American Heart Association (AHA 2017) table and provides a cue if your reading falls into one of the stages that could potentially indicate increased risk. See page 21 for more information on this feature.

The appearance of the  icon indicates that a pulse irregularity consistent with an irregular heartbeat was detected during measurement. Refer to page 22 for more information on the Irregular Heartbeat Detector.

NAME/FUNCTION OF EACH PART

Accessories for unit:

2 "AAA" 1.5V alkaline batteries and a carrying case are included with the unit.

DISPLAY EXPLANATIONS

Date/Time
Indicator

12/28/19
12:35 PM

Systolic
Pressure

115

Wrist
Positioning
Guide

Diastolic
Pressure

75

Heart Rate

72
Pul
/min

Irregular
Heartbeat
Detector

Display Symbols:


	User 1: Appears when the monitor is operated by User 1.
	User 2: Appears when the monitor is operated by User 2.
	Low Battery Symbol: Appears when batteries should be replaced, or +/- polarities have been positioned incorrectly.
	Pulse Symbol: Shows heart rate in beats per minute.
	Wrist Positioning Guide: See page 17 for more information.
	Irregular Heartbeat Detector: Irregular heartbeat detected during measurement. See page 22 for more information.
	Risk Category Index: Indicates the blood pressure level. See page 21 for more information.
AVG.3	Memory Average: Displays average of last 3 readings.
AVG.3 AM	AM Averaging: Indicates the reading being displayed is an average from the last 3 mornings measurements, within the last 7 days. AM is defined as 4:00 AM – 12:00 PM.
AVG.3 PM	PM Averaging: Indicates the reading being displayed is an average from the last 3 evening measurements, within the last 7 days. PM is defined as 6:00 PM – 2:00 AM.
	Excessive Body Motion Detector: Appears when excessive body movement, especially of the wrist the cuff is worn, is detected during the measurement. See page 23 for more information. NOTE: The measured blood pressure reading may not be accurate if the icon is displayed.

If any of the following letters and numbers appear in the area that systolic pressure should be displayed or the wrist icon appears, an error has occurred with your reading. See Troubleshooting section of this manual for more information.

E1	Incorrect cuff placement: Correct and measure again.
E2	Error Determining Measurement Data: Relax, do not move, and measure again.
E3	Pulse signal not detected: Measure again.
E4	Error Determining Measurement Data: Relax for 3 minutes, measure again.
EE	System Error: Turn off monitor and measure again. If EE still appears on the display followed by any other character (ex: EE4 , EEA , etc.), please contact Consumer Relations.
Out	Pressure Measurement Out of Range: Switch the unit off to clear, then measure again.

INSTALLING BATTERIES

1. Slide off the battery cover to open the battery compartment.
2. Install or replace the 2 “AAA” alkaline batteries according to the indications inside the battery compartment.

NOTE:

When removing batteries, pull black ribbon forward.

3. Replace the battery cover.

Replace the batteries if:

- The low battery symbol appears on the display.
- When any button is pressed and nothing is displayed on the screen.

NOTE:


- Date and time will need to be reset if batteries are removed or replaced.
- Replace all batteries at one time (as simultaneous set). Do not mix old and new batteries.
- Use only 1.5V “AAA” alkaline batteries. Do not mix alkaline, standard (carbon-zinc) or rechargeable (Ni-Cad, Ni-MH, etc.) batteries.
- When installing batteries, observe proper +/- polarities. Incorrect installation of batteries may cause damage to the unit.

- Remove batteries when unit is not in use for extended periods of time.
- When the batteries are removed, the measurement values stored in memory are retained. However, the date and time must be reset.
- Clean contacts on battery and in battery compartment with a soft dry cloth each time you install batteries.
- Keep batteries away from children as they could pose a choking hazard.
- Recycle or dispose of properly in accordance with all local, state, province, and country regulations.
- **Batteries are hazardous waste. DO NOT dispose of them together with the household garbage.**
- **DO NOT dispose of batteries in fire. Batteries may explode or leak.**


DATE & TIME SET PROCEDURE

It is important to set the Date & Time before using your blood pressure monitor, so that a time stamp can be assigned to each record that is stored in the memory.

1. To set the date and time, press the **“Date & Time Set”** ⌚ button.

2. The display will show a blinking number representing the MONTH. Change the MONTH by pressing the **“MEM”** button. Each press will increase the number by one in a cycling manner. Press the ⌚ button again to confirm the entry. The screen will then show a blinking number representing the DAY.
3. Change the DAY, YEAR, HOUR, & MINUTE as described in step 2 above, using the **“MEM”** to change the numbers and the ⌚ button to confirm the entries.

Repeat process to set DAY, YEAR, HOUR and MINUTE.

IMPORTANT: IF YOU ARE USING THIS BLOOD PRESSURE MONITOR FOR THE FIRST TIME, PLEASE REMOVE THE PROTECTIVE FILM FROM THE DISPLAY SCREEN.

APPLYING THE CUFF

1. Remove all watches, wrist jewelry, etc. prior to attaching the wrist monitor. Clothing sleeves should be rolled up and the cuff should be wrapped around bare skin for correct measurements.
2. Apply cuff to left wrist with palm facing up as shown in Fig. A.
3. Make sure the edge of the cuff is about $1/4''$ – $1/2''$ (1cm–1.5cm) from the palm as shown in Fig. B.
4. In order to ensure accurate measurements, fasten the hook and loop strap securely around your wrist so there is no extra space between the cuff and the wrist. If the cuff is not wrapped tight enough, the measurement values will not be accurate.
5. If your physician has diagnosed you with poor circulation in your left arm, place the cuff around your right wrist as shown in Fig. C.

Fig. A


Fig. B



Fig. C

NOTE:

- This device should not be used when your wrist is wounded or injured.
- If it is not possible to use the cuff on your left wrist, it can also be used on your right wrist. However, all measurements should be made using the same wrist.

CORRECT MEASUREMENT POSTURE

1. Place your elbow on a table so that the cuff is at the same level as your heart as shown in Fig. A. Relax your entire body, especially between your elbow and fingers.

NOTE:

Your heart is located slightly below your armpit.

2. If the cuff is not at the same level as your heart or if you cannot keep your arm completely still throughout the reading, use a soft object, such as a folded towel, to support your wrist as shown in Fig. B.
3. Turn your palm upwards.
4. Sit upright in a chair with feet flat on the floor and take 5-6 deep breaths. Avoid leaning back while the measurement is being taken as shown in Fig. C.

Fig. A


Fig. B


Fig. C

WRIST POSITIONING GUIDE

- When using wrist blood pressure monitors, it is important to take the measurement while relaxing your wrist at heart level. This will help ensure an accurate reading will be taken. Because it is difficult to find the ideal location for a wrist measurement, we have included an optional Wrist Positioning Sensor that will assist in directing your wrist to the ideal location for a blood pressure measurement.
- After you press “**START/STOP**” button, the display will illuminate with different icons that are designed to help you move your wrist. Once the ideal location is found, the wrist symbol  will flash twice and measurement will begin. Keep your arm still until measurement is completed.

SYMBOL	ACTION
	Move your wrist up.
	Move your wrist down.
	Rotate your wrist to your left.
	Rotate your wrist to your right.
	Wrist is in ideal position, do not move.
	Wrist was <u>NOT</u> in ideal position, this symbol will be saved in the measurement memory if this occurs.


- Even if the monitor is not in the ideal position after 10 seconds, the measurement will still begin. If this occurs, the wrist symbol  will appear and be stored in the memory with the measurement.

NOTE:


- Due to differences in individual sizes and physique, this feature may not be helpful in all cases and you may wish to turn this feature OFF. If you feel the suggested wrist position does not match your heart level, please turn this feature OFF and follow your judgment.

WRIST POSITIONING GUIDE SETTINGS

1. When the monitor is off, to turn the Wrist Positioning Guide OFF or ON, press and hold the **“START/STOP”** button for 3 seconds.

Wrist Positioning
Guide is ON

Wrist Positioning
Guide is OFF

NOTE:

- By default, the Wrist Positioning Guide is ON.
- When the Wrist Positioning Guide is OFF, the instruction will not appear at the start of the measurement.

MEASUREMENT PROCEDURE

IMPORTANT NOTES:

- Blood pressure changes with every heartbeat and is in constant fluctuation throughout the day.
- Blood pressure measurement can be affected by the position of the user, his or her physiologic condition, and other factors.

Before Measurement:


To help ensure an accurate reading, follow these instructions before taking a measurement:

- Wait 1 hour after exercising, bathing, eating, drinking beverages with alcohol or caffeine, or smoking to measure blood pressure.
- Sit quietly and rest for 15 minutes.
- Stress raises blood pressure. Avoid taking measurements during stressful times.
- Take your reading in a comfortable environment as measurements can be affected by hot or cold temperatures. Take your blood pressure at normal body temperature.


During Measurement:

- DO NOT talk or move your arm or hand muscles.
- DO NOT cross your legs. Sit with feet flat on the floor.
- DO NOT touch cuff or monitor during measurement.

1. Press the “1&2” button and select User 1 or User 2.
2. With the cuff wrapped around your wrist, press the “**START/STOP**” button. Do not inflate the cuff unless it is wrapped around your wrist. All digits will light up, checking the display functions. The checking procedure will be completed after about 3 seconds.

3. After all symbols disappear, the display will show “0.” The monitor is ready to measure and will automatically inflate the cuff to begin measurement.
4. When the measurement is completed, the cuff will deflate entirely. Systolic pressure, diastolic pressure, and pulse will be shown simultaneously on the LCD screen. **The measurement is then automatically stored into memory.**

NOTE:

- This monitor automatically turns off approximately **1 minute** after last operation. You may also press the “**START/STOP**” button to turn the unit off.
- To interrupt the measurement, you may press the “**START/STOP**” (recommended) or the “**MEM**” or “**Date & Time Set**” ⌚ buttons. The cuff will deflate immediately after this button is pressed.

RISK CATEGORY INDEX

This monitor comes equipped with a Risk Category Indicator that automatically compares each reading to defined levels established by the American Heart Association® (AHA 2017) as described earlier in this manual, and provides a helpful cue if your reading falls into one of the stages that could potentially indicate increased risk. Please note that cues provided by this monitor are only intended to assist you in using this table.


The table is only provided for your convenience to help you understand your noninvasive blood pressure readings as it relates to the AHA 2017 information. They are not a substitute for a medical examination by your physician. It is important for you to consult with your physician regularly. Your physician will tell you your normal blood pressure range as well as the point at which you may actually be considered to be at risk.

Blood Pressure Category	Systolic mmHg (upper number)		Diastolic mmHg (lower number)	Indicator Color
Normal	<120	and	<80	Green
Elevated	120–129	and	<80	Yellow
High Blood Pressure (hypertension) Stage 1	130–139	or	80–89	Red
High Blood Pressure (hypertension) Stage 2	140–180	or	90–120	
Hypertensive Crisis (consult your doctor immediately)	>180	and/ or	>120	

*Source: AHA 2017

IRREGULAR HEARTBEAT DETECTOR

The appearance of the  icon indicates that a pulse irregularity consistent with an irregular heartbeat was detected during measurement. Usually this is not a cause for concern. However, if the symbol appears often, we recommend you seek medical advice. Please note that the device does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage.

Irregular Heartbeat Detector

Movement, shaking, or talking during the measurement can result in pulse irregularities that may cause the appearance of this icon. Therefore, it is of great importance to not move or talk during measurement.

During each measurement, the blood pressure monitor will keep a record of all the pulse intervals and calculate the average value of them. If there are two or more pulse intervals with 25% or greater difference, or if there are four or more pulse intervals with 15% or greater difference from the average heartbeat interval, then the irregular heartbeat symbol  will appear on the display with the measurement result.


IMPORTANT INFORMATION:

This blood pressure monitor is not designed for use by people with arrhythmias nor for diagnosing or treating an arrhythmia problem. As a safeguard, we recommend that if you have arrhythmias such as atrial or ventricular premature beats and atrial fibrillation or any other special conditions, you should check with your physician before using your blood pressure monitor.

EXCESSIVE BODY MOTION DETECTOR

The “” icon is displayed if body movement is detected during measurement, especially the movement of the wrist the blood pressure monitor is worn on.

Excessive
Body Motion
Detector

NOTE:

The measured blood pressure reading may not be accurate if the icon is displayed.


RECALLING VALUES FROM MEMORY

This monitor can be used by two individuals. Each user can store up to 100 measurements.

This monitor features an advanced memory mode to provide you with a variety of options for reviewing your measurement history, including AM/PM averaging. Reviewing your morning (AM) and nighttime (PM) measurements can provide important information about how your blood pressure changes throughout the day.

To use the memory function:


1. Press the “” button and select User 1 or User 2.

2. Press the “**MEM**” button to access the memory.


3. The monitor will first display the calculated average applied to the last 3 memories.


NOTE:

- If the records are less than 3, it will display the latest record.
- The average will not include readings where excessive body motion (Ⓐ) occurred, or readings where the wrist is not in the ideal position (👉).

4. Press the **“MEM”** button again and the monitor will display the average of the last 3 AM readings within the last 7 days. Press the **“MEM”** button again to display the last 3 PM readings within the last 7 days.

AM Averaging


PM Averaging

NOTES:

- If there are fewer than three AM or PM measurements, this average will not be displayed.
 - AM is defined as 4:00 AM – 12:00 PM.
 - PM is defined as 6:00 PM – 2:00 AM.
5. Every new press of the **“MEM”** button will recall a previous reading. This latest reading will be recalled first.


NOTE:

If the wrist positioning guide is ON before the measurement, and the monitor is not in the ideal position, the 👉 icon will be displayed with the saved readings.

CLEARING VALUES FROM MEMORY


1. Press the “1 Δ /2 Δ ” button to select User 1 or User 2.

2. Press the “MEM” button to enter into memory recall mode.

3. Press and hold the “MEM” and the “Date & Time Set” \odot buttons at the same time and the data for the selected user will be erased automatically.

NOTE:

Once deleted, your readings cannot be restored.

IMPORTANT NOTES REGARDING YOUR BLOOD PRESSURE MEASUREMENT

- It is suggested that you take your measurements at the same time each day and use the same wrist for consistency.
- Users should wait a minimum of 5 minutes before taking additional measurements. More time may be necessary depending upon your physiology.
- The measurement results that users receive are for reference only. If users have any blood pressure concerns, please consult a physician.
- Once inflation reaches 300 mmHg, the unit will deflate automatically for safety reasons.
- This product is not suitable for people with arrhythmias.
- This device may have difficulty determining the proper blood pressure for pregnant women and for users with irregular heartbeat, diabetes, poor circulation of blood, kidney problems, or for users who have suffered a stroke.

CARE, MAINTENANCE, & CLEANING

- Clean the blood pressure monitor body and cuff carefully with a slightly damp, soft cloth. Do not press. Do not wash cuff or use chemical cleaner on it. Never use thinner, alcohol, or petrol (gasoline) as cleaner.
- Leaky batteries can damage the unit. Remove the batteries when the unit will not be used for a long time.
- Follow local ordinances and recycling instructions regarding disposal or recycling of the device and device components, including batteries.
- If the unit is stored near freezing, allow it to acclimate to room temperature before use.
- This blood pressure monitor is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems with this device, please contact HoMedics Consumer Relations (contact information on warranty page).
- DO NOT immerse the unit in water as this will result in damage to the unit.
- DO NOT subject the monitor or cuff to extreme temperatures, humidity, moisture, or direct sunlight. Protect from dust.
- DO NOT fold the cuff tightly.
- DO NOT disassemble the monitor or cuff. If in need of repair, refer to the Warranty section of this manual.
- DO NOT subject the monitor to extreme shock (do not drop on floor).

- DO NOT inflate the cuff unless wrapped around wrist.
- DO NOT wrap the cuff around body parts other than your wrist.
- DO NOT drop or insert any object into any opening.
- This monitor may not meet its performance specifications if stored or used outside of these temperature and humidity ranges:

Storage/Transportation Environment

Temperature:

-4°F~140°F (-20°C~60°C)

Humidity:

≤ 93% RH

Operating Environment

Temperature:

41°F~104°F (5°C~40°C)

Humidity:

15%~90% RH

Atmospheric Pressure:

700 hPa to 1060 hPa

POTENTIAL FOR ELECTROMAGNETIC INTERFERENCE

To avoid inaccurate results caused by electromagnetic interference between electrical and electronic equipment, do not use the device near a cell phone or microwave oven. For most wireless communication devices, it is recommended to maintain a distance of 10.8 feet (3.3m) in order to avoid electromagnetic interference. This device complies with Part 15 of the FCC Rules.

Operation is subject to the following 2 conditions: (1) This device may not cause harmful interference, and (2) This device must accept any interference received, including interference that may cause undesired operation.

FEDERAL COMMUNICATIONS COMMISSION COMPLIANCE STATEMENT

Changes or modifications to this equipment not expressly approved by the manufacturer could void the user's authority to operate the equipment.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

ELECTROMAGNETIC COMPATIBILITY (EMC)

• Guidance and manufacturer's declaration – electromagnetic emissions

The device is intended for use in the electromagnetic environments listed below, and should only be used in such environments:

Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	RF energy is used only to maintain device's operation. Therefore, its RF emissions are so low that it's not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	The device is suitable for use in all establishments, including domestic establishments, and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

• Guidance and manufacturer's declaration – electromagnetic immunity

The device is intended for use in the electromagnetic environments listed below, and should only be used in such environments:

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	power supply lines: ± 2 kV input/output lines: ± 1 kV	power supply lines: ± 2 kV	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	line(s) to line(s): ± 1 kV line(s) to earth: ± 2 kV 100 kHz repetition frequency	line(s) to line(s): ± 1 kV 100 kHz repetition frequency	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°: 0% UT; 1 cycle 70% UT; 25/30 cycles 0% UT; 250/300 cycles	0% UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°: 0% UT; 1 cycle 70% UT; 25 cycles 0% UT; 250 cycles	Mains power quality should be that of a typical commercial or hospital environment.
Power frequency (50 or 60 Hz) magnetic field IEC 61000-4-8	30 A/m 50 or 60 Hz	30 A/m 50 or 60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: UT is the a.c. mains voltage prior to application of the test level.

• **Guidance and manufacturer's declaration – electromagnetic immunity**

The device is intended for use in the electromagnetic environments listed below, and should only be used in such environments:

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance	
Conducted RF IEC 61000-4-6	150 kHz to 80 MHz 3 Vrms to 6 Vrms (in ISM and amateur radio bands) 80% Am at 1 kHz	150 kHz to 80 MHz 3 Vrms to 6 Vrms (in ISM and amateur radio bands) 80% Am at 1 kHz	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 0.35 \sqrt{P}$ $d = 1.2 \sqrt{P}$	
Radiated RF IEC 61000-4-3	10 Vrms, 80% Am at 1 kHz	10 Vrms, 80% Am at 1 kHz	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$ 800 MHz to 2.7 GHz $d = 2.3 \sqrt{P}$	where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer, and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range. b Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1

At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2

These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.

b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

• **Recommended separation distances between portable and mobile RF communication equipment and the device.**

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = 3.5 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.7 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

• **Guidance and manufacturer's declaration – electromagnetic immunity**

The device is intended for use in the electromagnetic environments listed below. The customer or the user of the device should only be used in such environments:

Immunity test	Test Frequency (MHz)	Band a) (MHz)	Service a)	Modulation b)	Modulation b) (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
Radiated RF IEC 61000-4-3 (Test specifications for ENCLOSE PORT IMMUNITY to RF wireless communications equipment) Radiated RF IEC 61000-4-3	385	380-390	TETRA 400	Pulse modulation b) 18 Hz	1.8	0.3	27
	450	430-470	GMRS 460, FRS 460	FM c) ± 5 kHz deviation 1 kHz sine	2	0.3	28
	710	704-787	LTE Band 13, 17	Pulse modulation b) 18 Hz	0.2	0.3	9
	745						
	780						
	810	800-960	GSM 800/900; TETRA 800; DECT; 820, CDMA 850, LTE B 5	Pulse modulation b) 18 Hz	2	0.3	28
	870						
	930						
	1720	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation b) 217 Hz	2	0.3	28
	1845						
	1970						
	2450	2400-2570	Bluetooth® WLAN, 802.11 B/G/N, RFID 2450, LTE Band 7	Pulse modulation b) 21 Hz	2	0.3	28
	5240	5100-5800	WLAN 802.11a/n	Pulse modulation b) 21 Hz	0.2	0.3	9
	5500						
5785							

NOTE: If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

a) For some services, only the uplink frequencies are included. b) The carrier shall be modulated using a 50% duty cycle square wave signal. c) As an alternative to FM modulation, 50% pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

The MANUFACTURER should consider reducing the minimum separation distance, based on RISK MANAGEMENT, and using higher IMMUNITY TEST LEVELS that are appropriate for the reduced minimum separation distance. Minimum separation distances for higher IMMUNITY TEST LEVELS shall be calculated using the following equation:

$$E = 6 \sqrt{\frac{P}{d}}$$

Where P is the maximum power in W, d is the minimum separation distance in m, and E is the IMMUNITY TEST LEVEL in V/m.

TROUBLESHOOTING

If any abnormality arises during use, please check the following points:

SYMPTOMS	POSSIBLE CAUSES	CORRECTION
Unit does not turn on when the "START/STOP" button is pushed.	Batteries have run down.	Replace them with 2 new "AAA" alkaline batteries.
	Battery polarities have been positioned incorrectly.	Re-insert the batteries in the correct positions.
[ERROR] E1 error symbol shown on display.	Cuff has been placed incorrectly or there is an air circuit abnormality.	Wrap the cuff properly so that it is positioned correctly. Measure again.
[ERROR] E2 error symbol shown on display.	Did you talk or move during measurement?	Keep wrist steady during measurement. Measure again.
	Shaking of the wrist with the cuff on.	
[ERROR] E3 error symbol shown on display.	Pulse signal not detected.	Measure again.
[ERROR] E4 error symbol shown on display.	Error determining measurement data.	Relax for 3 minutes, measure again.
[ERROR] EE error symbol shown on display.	System error.	Turn off monitor and measure again. If EE still appears on the display followed by any other character (ex: EE4, EEA, etc.), please contact Consumer Relations.
out error symbol shown on the display	Pressure measurement out of range.	Switch the unit off to clear, then measure again.

 <p>Excessive body motion detector symbol shown on display.</p>	<p>Excessive body movement, especially of the wrist the cuff is worn, occurred during the measurement. Ex: talking, moving, or shaking of the arm with the cuff on while measurement.</p>	<p>Measure again. Keep arm steady during measurement.</p>
	<p>Cuff may not be wrapped correctly.</p>	<p>Ensure wrist cuff is wrapped correctly. See "Applying the Cuff" on page 15 for more information. Please consult your health care professional if you have questions on how to properly wear the wrist cuff.</p>

Note: If the unit still does not work, contact HoMedics Consumer Relations. Under no circumstance should you disassemble or attempt to repair the unit by yourself. Contact information for HoMedics Consumer Relations Department can be found on the warranty page.

SPECIFICATIONS

Power Source:	Two 1.5V DC “AAA” Alkaline batteries
Measurement Method:	Oscillometric
Measurement Range:	Rated cuff pressure: 0 mmHg~299 mmHg (0kPa~39.9 kPa) Measurement pressure: SYS: 60 mmHg~230 mmHg (8.0 kPa~30.7 kPa) DIA: 40 mmHg~130 mmHg (5.3 kPa~17.3 kPa) Pulse value: 40~199 beats/minute
Accuracy:	Pressure: ± 3 mmHg Pulse: $\pm 5\%$ of reading
Pressure Sensor:	Semi-conductor
Inflation:	Pump driven
Deflation:	Automatic Pressure release valve
Memory Capacity:	200 memories
Auto Shutoff:	1 minute after last button operation
Battery Life:	Approx. 150 measurements
Operation Environment:	Temperature: 41°F~104°F (5°C~40°C) Humidity: 15% ~ 90% RH Atmosphere pressure range: 700 hPa to 1060 hPa
Storage/Transportation Environment:	Temperature: -4°F~140°F (-20°C~60°C) Humidity: $\leq 93\%$
Weight:	Approx. 0.26 lbs (115 g) (without batteries)

Wrist Circumference:	5.3"-8.4" (13.5cm-21.5cm)
Accessories:	(2) "AAA" alkaline batteries, instruction manual, carrying case.
Dimensions:	2.75"(L) x 3.58"(W) x .93"(H) 70mm(L) x 91mm(W) x 23.5mm(H)
Degree of Protection	Type BF applied part
Device Classification	Internally Powered ME Equipment
IP Classification	IP22: The first number 2: Protected against solid foreign objects of 12.5mm Ø and greater. The second number: Protected against vertically falling water drops when enclosure tilted up to 15°. Vertically falling drops shall have no harmful effects when the enclosure is tilted at any angle up to 15° on either side of the vertical.
Software Version	A01
	Follow instructions for use.

Made exclusively for

Walgreens

**For service or repair,
do not return this unit
to the retailer. Contact
HoMedics Consumer
Relations at:**

**Email:
cservice@
homedics.com**

**Phone:
1-800-466-3342
Business Hours:
8:30am-7pm ET
Monday-Friday**

**Distributed by
HOMEDICS®**

**HoMedics USA, LLC
3000 N Pontiac Trail
Commerce Township, MI
48390**

Printed in China

LIMITED 5-YEAR WARRANTY

HoMedics sells its products with the intent that they are free of defects in manufacture and workmanship for a period of 5 years from the date of original purchase, except as noted below. HoMedics warrants that its products will be free of defects in material and workmanship under normal use and service. This blood pressure monitor meets the simulated measurement cycles test requirement per EN1060-3, part 8.10. This warranty extends only to consumers and does not extend to Retailers.

To obtain warranty service on your HoMedics product, contact a Consumer Relations Representative by telephone at 1-800-466-3342 for assistance. Please make sure to have the model number of the product available.

HoMedics does not authorize anyone, including, but not limited to, Retailers, the subsequent consumer purchaser of the product from a Retailer or remote purchasers, to obligate HoMedics in any way beyond the terms set forth herein. This warranty does not cover damage caused by misuse or abuse; accident; the attachment of any unauthorized accessory; alteration to the product; improper installation; unauthorized repairs or modifications; improper use of electrical/power supply; loss of power; dropped product; malfunction or damage of an operating part from failure to provide manufacturer's recommended maintenance; transportation damage; theft; neglect; vandalism; or environmental conditions; loss of use during the period the product is at a repair facility or otherwise awaiting parts or repair; or any other conditions whatsoever that are beyond the control of HoMedics. This warranty is effective only if the product is purchased and operated in the country in which the product is purchased. A product that requires modifications or adoption to enable it to operate in any other country than the country for which it was designed, manufactured, approved and/or authorized, or repair of products damaged by these modifications is not covered under this warranty.

THE WARRANTY PROVIDED HEREIN SHALL BE THE SOLE AND EXCLUSIVE WARRANTY. THERE SHALL BE NO OTHER WARRANTIES EXPRESS OR IMPLIED INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS OR ANY OTHER OBLIGATION ON THE PART OF THE COMPANY WITH RESPECT TO PRODUCTS COVERED BY THIS WARRANTY. HOMEDICS SHALL HAVE NO LIABILITY FOR ANY INCIDENTAL, CONSEQUENTIAL, OR SPECIAL DAMAGES. IN NO EVENT SHALL THIS WARRANTY REQUIRE MORE THAN THE REPAIR OR REPLACEMENT OF ANY PART OR PARTS WHICH ARE FOUND TO BE DEFECTIVE WITHIN THE EFFECTIVE PERIOD OF THE WARRANTY. NO REFUNDS WILL BE GIVEN. IF REPLACEMENT PARTS FOR DEFECTIVE MATERIALS ARE NOT AVAILABLE, HOMEDICS RESERVES THE RIGHT TO MAKE PRODUCT SUBSTITUTIONS IN LIEU OF REPAIR OR REPLACEMENT.

This warranty does not extend to the purchase of opened, used, repaired, repackaged and/or resealed products, including but not limited to sale of such products on internet auction sites and/or sales of such products by surplus or bulk resellers. Any and all warranties or guarantees shall immediately cease and terminate as to any products or parts thereof which are repaired, replaced, altered, or modified, without the prior express and written consent of HoMedics.

This warranty provides you with specific legal rights. You may have additional rights which may vary from state to state. Because of individual state regulations, some of the above limitations and exclusions may not apply to you.

**©2019 HoMedics USA, LLC. All rights reserved.
HoMedics is a registered trademark of HoMedics, LLC.
GentleInflate is a trademark of HoMedics, LLC.**

American Heart Association is a registered trademark of the American Heart Association.

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©2019 HoMedics USA, LLC. Todos los derechos reservados.
HoMedics es una marca registrada de HoMedics, LLC.
GentleInflate es una marca comercial de HoMedics, LLC.